



CPAL

Central Pennsylvania Alliance Laboratory

Technical Bulletin Number 20

Prealbumin Reference Range Change

In the month of November, CPAL will be converting Prealbumins to a new testing platform. With this change, there will be a modification in the reference range for Prealbumins.

The current reference range is 15.0 to 32.0 mg/dL.

The new reference range will be 18.0 to 35.7 mg/dL.

When this change goes into effect, a comment will be attached to the results which states:

Note reference range change effective (date).

Issued on: January 12, 2004

For question about this, and other, information, call Central Pennsylvania Laboratory at 1-888-480-1422